

Technical Note: Using Ethernet Switch for CCU

Introduction

In general, a long cable is required to connect the common control unit (CCU) to the local network or PC. The long cable may reduce the signal significantly. An Ethernet Switch serves as signal amplifier, and also provides more parallel connections to several devices.

Connect the CCU via Ethernet Switch to the network

Ethernet Switch delivered with the shipment: 5-Port 10/100 Fast Ethernet Switch FS 605.

1. Plug one end of an Ethernet cable to local network and the other end to one of the Ethernet port on the switch.
2. For each device (CCU, PC, emXX, etc., up to 4 devices), plug one end of the Ethernet cable to the device, and the other end to one of the Ethernet port on the switch.
3. Connect the power adapter to the switch, and then plug the adapter to the power supply.

For more information the Ethernet Switch, please refer to the switch Installation Guide.

In case of any further questions, please contact us under:

IBA Dosimetry GmbH
Service Department
Bahnhofstrasse 5
DE-90592 Schwarzenbruck, Germany
phone: +49 9128 607 38
fax: + 49 9128 607 10
mailto: service@iba-group.com