
## **X-LITE**

## QA Tool for Light Field Alignment Checks



Light up your radiation fields.

#### **Fast and Precise Alignment**

Correct alignment of the light field with the radiation field is essential for reliable treatment set-up. According to national and international recommendations, e.g. NACP (Acta Radiologica Oncology 19 Fasc 1, 1980) and IAEA (Technical Report Series No 277, 1987), this alignment should be checked at least once a week and should correspond to within  $\pm$  2 mm at SSD 100 cm.

The X-LITE with a fluorescent plate helps you to align the accelerator radiation fields on the treatment table quickly and directly. The plate is activated by ionizing radiation, producing a green fluorescence which is clearly visible for a few minutes after the radiation is switched off. The  $5 \times 5$ ,  $10 \times 10$ ,  $15 \times 15$  and  $20 \times 20$  cm² active areas are permanently marked with field scales in centimetres.

A Protective Red Filter prevents activation of the phosphor by the light field from the gantry head or by ambient light. Irradiation photons and electrons penetrate the filter to activate the plate.

Using X-LITE is Easy, and alignment checks can be part of the daily routine. With the filter in place, align the light field according to the scale on X-LITE. Put any required build-up material on the plate and irradiate. Remove the filter from X-LITE and check the radiation field alignment.

It is not necessary to wait for the afterglow to fade between irradiations, since the contrast is high as long as the following irradiation uses the same dose or higher.



# **Technical Specifications**

X-LITE	
Active area:	23 cm x 23 cm
Field scale:	5 cm x 5 cm, 10 cm x 10 cm, 15 cm x 15 cm
	and 20 cm x 20 cm
Deviation scale:	$\pm$ 5 mm in steps of 1 mm from each field scale
Scale accuracy:	± 1 mm
Maximum dose:	6 Gy/exposure
Daylight filter:	red, removable
Material:	plexiglass / polycarbonate
Operating temperature:	15 – 45°C
Dimensions:	320 mm (L) x 276 mm (W) x 14 mm (H)
Weight:	1.4 kg

Technical data is subject to change without prior notice.

### Light up your radiation fields.

Contact details:

dosimetry-info@iba-group.com

**Europe, Middle East, Africa** 

IBA Dosimetry GmbH Bahnhofstr. 5

90592 Schwarzenbruck, Germany

Tel.: +49 9128 607 0 Fax: +49 9128 607 10

North America, Latin America

IBA Dosimetry America 3150 Stage Post Drive, Suite 110 Bartlett, TN 38133, USA

Tel.: +1 901 386 2242 Fax: +1 901 382 9453

**Asia Pacific** 

IBA Dosimetry Asia Pacific No.6, Xing Guang Er Jie Beijing OPTO-mechatronics Industrial Park (OIP), Tongzhou District Beijing 101111, China

Tel.: +86 10 8080 9288 Fax: +86 10 8080 9299

