

Compression-pressure meter for mammography exposure units

General explanation / description:

The compression-pressure meter for mammography exposure units is made for the measurement of the constancy and accuracy of the compression of mammography units how described in the "European Guidelines for quality assurance ..." for the yearly test.

The advantages of this measure device are in the easy operation and the very exact measure results over the whole measure plate.

The big LC-Display give the possibility of quickly control of the measured results.

Technical data

- Measure range: 200 g to 25 kg
- Resolution of the LCD: 10 g
- Power supply: AC/DC adapter
12 V 1.5 A
- Temperature range: 15°C – 35°C
- Admission: OIML 4000d for
Scales class III
and class IIII



Status: September 2008