


## Quality Assurance in Radiodiagnostics

### Filter Sets for HVL measurements |

Aluminum attenuator set for **HVL** measurements  
for X-ray units working in the range of 40 - 150 kV



<b>Dimensions:</b>	100 x 100 mm each
<b>Purity of Al:</b>	$\geq 99.5 \%$
<b>Inaccuracy of thickness:</b>	$\pm 10 \%$
<b>Set consisting of:</b>	5 filter plates of 0.1 mm Al 2 filter plates of 0.5 mm Al 5 filter plates of 1.0 mm Al 2 filter plates of 2.0 mm Al

Aluminum attenuator set for **HVL** measurements  
for for mammography units



<b>Dimensions:</b>	100 x 100 mm each
<b>Purity of Al:</b>	$\geq 99.99 \%$
<b>Inaccuracy of thickness:</b>	$\pm 10 \%$
<b>Set consisting of:</b>	7 filter plates of 0.1 mm Al