

## MiniTRACE C-10 Contamination Meter



- **Display unit:** Counts per second (cps)
- **Display- and measurement range:** 0.0 – 9990.0 cps
- **Sensitivity:** 0.15 cps/Bq (Co-60)
- **Display:** 6-digit LC-display, additional 5-digit alphanumeric display for system quality and alarm features
- **Alarm:** Acoustical and visual alarm
- **Buttons:** 2 buttons to switch on/off
- **Detector:** Geiger-Mueller pancake, active counter tube area 15.55 cm<sup>2</sup>, active diameter 44.5 mm, max. diameter 53.6 mm, entry window 2.0 mg/cm<sup>2</sup>, HV range 500 V
- **Grid:** 0.8 mm stainless steel, easily demountable, 80 % transparency
- **Energy supply/autonomy:** 2 Mignon batteries (type: LR6, AA, MN 1500), protected against wrong poling, battery life-time ca 2000 hrs

Status: July 2008