

## Technical Note FD19\_TN001\_081222

OmniPro Accept and OmniPro I'mRT

Module: Vidar film scanner option

### Vidar Scanners and Toolkit Overview

Overview:

Scanner	FW	Toolkit	Windows	ASPI/STI	Notes	
VXR-12	3.x 4.x	Filedate 1999-10-06	98	ASPI	Upgrade FW to 5.05	
			NT			
	-	2000	-	Not validated		
	-	XP	-			
VXR-12plus	5.x	Filedate 1999-10-06	98	ASPI	Upgrade FW to 5.05	
			NT	ASPI		
	-	2000	-	Not validated		
	-	XP	-	Not validated		
VXR-16	7.x- 12.x	Filedate 1999-10-06	98	ASPI	Upgrade FW to 13.0	
			NT	ASPI		
			-	2000		-
			-	XP		-
	13.0	3.1.9.3	98	ASPI		
			NT	ASPI		
			2000	STI		
			XP	STI		
VXR-16 Dosimetry Pro	7.x	Filedate 1999-10-06	98	ASPI	Upgrade FW to 13.0	
			NT			
			2000	-		
			XP	-		
	13.0	3.1.9.3	98	ASPI		
			NT	ASPI		
			2000	STI		
			XP	STI		
Dosimetry Pro Advantage	40.7	3.1.9.3	98	ASPI	Upgrade FW to 45.5	
			NT	ASPI		
			2000	STI		
			XP	STI		
	45.5	4.1.81.0	98	ASPI	Serial no 300485-	
			NT	ASPI		
			2000	STI		
			XP	STI		
	46.3	4.1.81.0	2000	STI	Serial no 310000-	
			XP	STI		



In case of any further questions, please contact us under:

IBA Dosimetry GmbH  
Service Department  
Bahnhofstrasse 5  
DE-90592 Schwarzenbruck  
Germany

Service hotline:  
phone: +49 9128 607 38  
fax: + 49 9128 607 10  
mailto: [service@iba-group.com](mailto:service@iba-group.com)