


VISICOIL™

LINEAR FIDUCIAL MARKER

new!



Designed for Interventional Radiology Approved by Interventional Radiology

VISICOIL's 0.50mm and 0.35mm markers are now available pre-loaded in 21 and 22 gauge 8cm and 15cm needles. Ideal for Interventional Radiology allowing for easier access with improved patient comfort and safety due to the smaller needle size.

IMPLANTATION SAFETY

- Reduced risk of pneumothorax
- Reduced risk of bleeding / edema
- Improved safety and comfort for patient

IMPLANTATION EFFICIENCY

- Shorter 8cm and 15cm needles ideal for standard bore CT
- Steerable needle aids in accurate marker placement
- Centimeter marks for depth placement



www.visicoil.com

Protect,
enhance
and save
lives

Iba

Safety, Efficiency and Clinical Confidence in Tumor Targeting

PATIENT SAFETY

- Conventional gold markers require a 17g needle
- VISICOIL IR needles are up to 50% smaller




- 17g
- 21g
- 22g

VISICOIL Needle Size

TREATMENT EFFICIENCY

- One marker - multiple reference points
- Flexible marker - visualize tissue deformation & motion



Pancreas: 2 VISICOIL 0.50mm x 1cm

CLINICAL CONFIDENCE

- Minimized artifact for accurate tumor targeting
- Non-migrating design



Liver: 1 VISICOIL 0.35mm x 1cm

NEW VISICOIL IR PRODUCT RANGE

VISICOIL SIZE	NEEDLE SIZE	PART NUMBER
0.50mm x 0.5cm	21 gauge x 8cm	SN-050-005-PL08
0.50mm x 0.5cm	21 gauge x 15cm	SN-050-005-PL15
0.50mm x 1.0cm	21 gauge x 8cm	SN-050-005-PL08
0.50mm x 1.0cm	21 gauge x 15cm	SN-050-010-PL15
0.35mm x 0.5cm	22 gauge x 8cm	SN-035-005-PL08
0.35mm x 0.5cm	22 gauge x 15cm	SN-035-005-PL15
0.35mm x 1.0cm	22 gauge x 8cm	SN-035-010-PL08
0.35mm x 1.0cm	22 gauge x 15cm	SN-035-010-PL15

VISICOIL is also available loose in various sizes and lengths to be used in the needle of your choice.

CONTACT DETAILS

NORTH AMERICA | LATIN AMERICA

RadioMed Corporation
3150 Stage Post Drive, Suite 110
Bartlett, TN 38133
Tel: +1 901 432 7202
Fax: +1 901 432 7206
Toll Free: 866 429 0922
visicoil@iba-group.com

EUROPE | MIDDLE EAST | ASIA PACIFIC

IBA Dosimetry GmbH
Bahnhofstr.5
90592 Schwarzenbruck, Germany
Tel: +49 9128 607 0
Fax: +49 9128 607 10
email: info@iba-dosimetry.com

PLEASE SEE BELOW CODE FOR IMPLANTATION INSTRUCTIONS



CALL CUSTOMER SERVICE FOR YOUR LOCAL REPRESENTATIVE
800 429 0922